GALEAS™/ Bladder

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Informed Genomics Limited

Laboratory No. 22262

**is accredited in accordance with International Standard ISO 15189:2012 -
Medical laboratories - Requirements for quality and competence**

This accreditation demonstrates technical competence for a defined scope specified in the schedule to this certificate, and the operation of a management system generally in accordance with the principles of ISO 9001 (refer to the Joint IAF-ILAC-ISO Communiqué on the Management system requirements of ISO 15189). The schedule to this certificate is an essential accreditation document and from time to time may be revised and reissued.

The most recent issue of the schedule of accreditation, which bears the same accreditation number as this certificate, is available from www.ukas.com.

This accreditation is subject to continuing conformity with United Kingdom Accreditation Service requirements

Matt Gantley, Chief Executive Officer
United Kingdom Accreditation Service

Initial Accreditation: 4 August 2021
Certificate Issued: 13 April 2023

Scan QR Code to
verify

UKAS is appointed as the sole national accreditation body for the UK by The Accreditation Regulations 2009 (SI No 3155/2009) and operates under a Memorandum of Understanding (MoU) with the Department for Business, Energy and Industrial Strategy (BEIS).

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
The logo for GALEAS™/BLADDER is displayed in white, uppercase letters on a dark blue rectangular background with a vertical gradient from dark blue on the left to teal on the right.

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Circulating Tumour DNA Detection By The Urine-Informed Analysis Of Archival Serum Samples From Muscle-Invasive Bladder Cancer Patients [Link para o Pubmed](#)

- Zhu, C. Z., Ting, H. N., Ng, K. H., & Ong, T. A. (2019). A review on the accuracy of bladder cancer detection methods. *Journal of cancer*, 10(17), 4038. [Link para o artigo](#)

- Alimi, Q., Hascoet, J., Manunta, A., Kammerer-Jacquet, S. F., Verhoest, G., Brochard, C., ... & Peyronnet, B. (2018). Reliability of urinary cytology and cystoscopy for the screening and diagnosis of bladder cancer in patients with neurogenic bladder: a systematic review. *Neurourology and Urodynamics*, 37(3), 916-925. [Link para o artigo](#)

- Ward DG, Gordon NS, Boucher RH, et al. Targeted deep sequencing of urothelial bladder cancers and associated urinary DNA: a 23-gene panel with utility for non-invasive diagnosis and risk stratification. *BJU Int*. 2019;124(3):532-544.

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
	Sensitivity	Specificity	PPV	NPV
pTa	87%	85%	75%	93%
T1	95%	85%	66%	98%
T2+	86%	85%	61%	96%
G1	78%	85%	47%	96%
G2	91%	85%	67%	97%
G3	91%	85%	78%	94%
NMIBC	90%	85%	83%	91%
MIBC	86%	85%	61%	96%

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
GALEAS/BLADDER

GALEAS Bladder Cancer Test

Patient ID: p11-83CA
Patient Name: Anon
Date of Birth: 01/01/1950
Sample ID: 006-090023
Received: 09/02/2023
Processed: 08/02/2023
Report Date: 05/02/2023 08:55
Sample Type: Urine Pellet
Clinician: Dr Dhye Deepen
Address: Coombs Ford Surgery
 Upper Downtown
 F14 8E8
 United Kingdom
Customer: Testing Laboratory
Address: Unit 5
 Quinton Business Park
 Quinton, Birmingham
 B32 7AF

Result Summary: Positive

A positive test result indicates that, at the time of GALEAS® Bladder testing, there is a high likelihood that cancer is present. A GALEAS® Bladder positive test result, at the time of testing, is determined by the presence of bladder cancer associated genomic variants (listed below) within the urine sample. The cancer association of the genomic markers used to inform the GALEAS® Bladder test have been reported in over 900 patient samples.

Variant Details

Gene	Ref/Alt	VAR
ARSA	c.175G>T	0.70%
SPRR1	c.265G>T	0.80%
CTHNA8	c.191G>T	0.49%
CTHNA8	c.191C>A	0.52%
PKC2A	c.248G>T	0.53%
PKC2A	c.353G>A	0.90%
PKC2A	c.302G>T	0.89%
PTPRF	c.398G>T	0.54%
PTPRF	c.156G>T	1.74%
PTPRF	c.394G>A	76.35%
YF11	PKCNA/G8	63.40%
CDKN2A	c.282G>T	0.70%
ESR1	c.126C>A	0.90%
ESR1	c.182C>A	1.89%
ESR1	c.270G>T	0.89%
TP53	c.335G>T	0.89%
TP53	c.332G>T	0.32%
TP53	c.876G>T	0.90%
TP53	c.574G>T	0.89%
TP53	c.838G>T	1.80%
TP53	c.567G>T	1.42%
TP53	c.548C>G	1.58%
TP53	c.539G>T	1.79%
TP53	c.827G>A	0.29%
TP53	c.323C>A	0.49%
TP53	c.326G>A	0.43%
TP53	c.840G>T	0.82%
TP53	c.150G>T	1.66%
PROX1	c.128G>T	0.74%

Test Description

GALEAS® Bladder is a next-generation sequencing-based test that assesses the detection and classification of bladder cancer associated genomic variants using urinary DNA specimens. The results of the test are used as an aid in identifying patients who may need further investigation due to the potential presence of bladder cancer.

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